

PubMed



Display Settings: Abstract



J Vis Exp. 2012 Nov 23;(69). pii: 4356. doi: 10.3791/4356.

Genome-wide gene deletions in *Streptococcus sanguinis* by high throughput PCR.

Ge X, Xu P.

The Philips Institute of Oral and Craniofacial Molecular Biology, **Virginia Commonwealth University**, Richmond, VA, USA.

Abstract

Transposon mutagenesis and single-gene deletion are two methods applied in genome-wide gene knockout in bacteria (1,2). Although transposon mutagenesis is less time consuming, less costly, and does not require completed genome information, there are two weaknesses in this method: (1) the possibility of a disparate mutants in the mixed mutant library that counter-selects mutants with decreased competition; and (2) the possibility of partial gene inactivation whereby genes do not entirely lose their function following the insertion of a transposon. Single-gene deletion analysis may compensate for the drawbacks associated with transposon mutagenesis. To improve the efficiency of genome-wide single gene deletion, we attempt to establish a high-throughput technique for genome-wide single gene deletion using *Streptococcus sanguinis* as a model organism. Each gene deletion construct in *S. sanguinis* genome is designed to comprise 1-kb upstream of the targeted gene, the *aphA-3* gene, encoding kanamycin resistance protein, and 1-kb downstream of the targeted gene. Three sets of primers F1/R1, F2/R2, and F3/R3, respectively, are designed and synthesized in a 96-well plate format for PCR-amplifications of those three components of each deletion construct. Primers R1 and F3 contain 25-bp sequences that are complementary to regions of the *aphA-3* gene at their 5' end. A large scale PCR amplification of the *aphA-3* gene is performed once for creating all single-gene deletion constructs. The promoter of *aphA-3* gene is initially excluded to minimize the potential polar effect of kanamycin cassette. To create the gene deletion constructs, high-throughput PCR amplification and purification are performed in a 96-well plate format. A linear recombinant PCR amplicon for each gene deletion will be made up through four PCR reactions using high-fidelity DNA polymerase. The initial exponential growth phase of *S. sanguinis* cultured in Todd Hewitt broth supplemented with 2.5% inactivated horse serum is used to increase competence for the transformation of PCR-recombinant constructs. Under this condition, up to 20% of *S. sanguinis* cells can be transformed using ~50 ng of DNA. Based on this approach, 2,048 mutants with single-gene deletion were ultimately obtained from the 2,270 genes in *S. sanguinis* excluding four gene ORFs contained entirely within other ORFs in *S. sanguinis* SK36 and 218 potential essential genes. The technique on creating gene deletion constructs is high throughput and could be easy to use in genome-wide single gene deletions for any transformable bacteria.

PMID: 23207952 [PubMed - indexed for MEDLINE]

Publication Types, MeSH Terms, Grant Support

LinkOut - more resources